

Complete Filter Cartridge Analyzer

Applications


Fully assembled filter cartridges are widely used in many industries including biotechnology, pharmaceutical, chemical, beverage, and food. Filtration efficiency in all these applications is determined by the pore structure of the complete filter cartridge. The Filter Cartridge Analyzer measures the bubble point, the mean flow pore diameter, and the pore distribution of the complete cartridge rather than a small sample of the filter media. The tester also measures the gas permeability.



The Filter Cartridge Analyzer

Principle

A wetting liquid is allowed to spontaneously fill the pores of the cartridge and air pressure is increased to empty pores and permit gas flow. Measured differential gas pressure and flow rates through a cartridge in wet and dry conditions yield various pore structure characteristics.



Principle of Cartridge Analyzer

The Instrument

The unique instrument design easily deals with very high volume of gas flowing through large cartridges and reduces pressure drop by eliminating narrow ducts, bends and constrictions. A tank for storage of gas under pressure sufficient for the test is supplied as part of the test equipment so that standard laboratory air supply is adequate for test execution. The sample chamber holds the cartridge between a fixed head and an adjustable head. By adjusting the position of the adjustable head cartridge of any length can be accommodated. A pneumatically operated piston applies sufficient pressure to seal the edges of the cartridge.



Sample Chamber for the Complete Cartridge

Capability

Pore diameter: The measured differential gas pressure yields the through pore throat diameter.

$$D = 4\gamma \cos \theta / p$$

D = pore diameter

γ = surface tension of wetting liquid

θ = contact angle of the liquid

p = differential gas pressure.

Bubble point: Computed from pressure for initiation of flow through wet sample.


Mean flow pore diameter: Computed from mean flow pressure at which wet curve and half-dry curves meet.

Pore distribution: Given in terms of distribution function, f.

$$f = -d[(F_w/F_d) \times 100] / dD$$

where F_w and F_d are wet and dry flow respectively.


Gas permeability: Computed from the gas flow rates using Darcy's law.



High gas flow rates through a filter cartridge measured as a function of differential pressure

Features

- ◆ Completely automated
- ◆ Windows based software, simple operation, and minimal operator involvement
- ◆ Only a few minutes for test execution
- ◆ Sintered metal, woven metal, polymeric, and ceramic cartridges can be tested
- ◆ Adequate safety precautions
- ◆ Normal laboratory air supply adequate to run the test



Pore distribution in a complete cartridge

Other Products

- | | |
|---|------------------------------|
| ◆ Porometers | ◆ Gas and Vapor Permeameters |
| ◆ Liquid Extrusion Porosimeters | ◆ Liquid Permeameters |
| ◆ Mercury/Nonmercury and Water Intrusion Porosimeters | ◆ Pycnometers |
| ◆ Gas Adsorption Sorptometers | ◆ Contract Testing Services |
| | ◆ Consulting Services |

Sale Rent Lease

Porous Materials, Inc.



20 Dutch Mill Road, Ithaca, New York 14850 USA

Toll Free US & Canada: 1-800-TALK-PMI Phone: (607)257-5544

Fax: (607)257-5639 Email: info@pmiapp.com Website: www.pmiapp.com

PMI Europe

Koningin Fabiolapark 45, BE 9820 Merelbeke, Belgium

Phone: +32 477 79 6011 Fax: +32 9 330 8544 Email: patrice.hellebaut@pmiappeurope.com